



MOLECULAR. IN MINUTES.™

# ID NOW™ STREP A 2


Now the fastest Strep A test  
without the need for culture  
confirmation.

- The fastest Strep A test – positive results in as few as 2 minutes<sup>1</sup>
- Enhanced molecular sensitivity eliminates the need for culture confirmation of negative results<sup>2</sup>
- Nearly 90% fewer false negatives than RADTs<sup>†,2,3</sup>
- CLIA waived<sup>2</sup>




# ID NOW™ Strep A 2

## SPEED — THE FASTEST STREP A TEST FOR EXPEDITING WORKFLOW



- Positive results detected in as few as 2 minutes, with 99% detected within 3 minutes<sup>1</sup>
- Negative results in 6 minutes<sup>2</sup>
- No incubation time or thermocycling required<sup>2</sup>

## PERFORMANCE — SUPERIOR SENSITIVITY ENHANCES CLINICAL CONFIDENCE



- Eliminates need for culture confirmation of negative results<sup>2</sup>
- Reduces unnecessary antibiotic usage

## EASE OF USE — SIMPLIFYING THE TEST PROCESS

- Room temperature storage
- CLIA waived for throat swab
- Integrated pipette eliminates manual pipetting
- Standardized respiratory trio (Influenza A & B, RSV, Strep A) across all sites of care
- Bi-directional connectivity to industry-leading AegisPOC™ platform facilitates satellite testing and remote management

## ORDER INFORMATION

PRODUCT NAME	PRODUCT CODE
ID NOW™ Strep A 2 – 24 Test Kit	734-000
ID NOW™ Strep A 2 – Control Swab Kit	734-080

**FOR MORE INFORMATION, CONTACT YOUR LOCAL ABBOTT REPRESENTATIVE OR VISIT [ABBOTT.COM/POCT](http://ABBOTT.COM/POCT).**

1. ID NOW™ Strep A 2 clinical trial data, held on file

2. ID NOW™ Strep A 2 Product Insert

3. Cohen JF, et al. Rapid antigen detection test for group A streptococcus in children with pharyngitis (Review). Cochrane Database of Systematic Reviews. 2016, Issue 7

†Based on the reported sensitivities of RADTs and ID NOW Strep A 2.

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. 120004442-03 09/19

